

ORTHODONTICS
KIEFERORTHOPÄDIE

NEWZED



SMAHA Needle holder
delicate tips

SMAHA Nadelhalter
feine Arbeitsenden

[140mm
175/20

slim and delicate
working tips – easier
access especially in
narrow, limited situations

Grazile Arbeitsspitzen –
ideal bei geringem
Bracketabstand



Distal End Cutter with TC inserts

Distalschneider
mit Hartmetalleinsätzen (TC)

⊗ Max. Ø 0.4mm / .016"

[150mm
632/16 TC

□ Max. .021x.025"

○ Ø .016"

[150mm
632/17 TC

□ Max. .021x.025"

○ Ø .016"

[130mm
632/18 TC

Rounded cutter head to ensure
pleasant treatment for the patient
and easier access for the user
to difficult to reach areas

Abgerundetes Arbeitsende,
angenehm für den Patienten und
ideal für den Anwender im
Bereich der hinteren Molaren



FON: +49 (0) 7463-7232

FAX: +49 (0) 7463-222

info@leibinger-dental.de

www.leibinger-dental.de